

**Ex Demo & Used Balance & Scales List.**  
**If you can't find what you need please contact us.**

**Precision Balance 300 g x 0.001 g**

Ex demo price £250 AML new price £355



Accuracy with economy from this new precision balance. A solidly built machine, it is able to weigh in multiple weighing units including Gram, Kilogram, Pound, Ounce, Carat, Tola (India), Troy Ounce, Pennyweight, Grain, Tael (3 units) & Centigram. One 200 g. and one 100 g. F2 grade calibration weights are included. Powered by supplied mains adaptor or 8 x AA size batteries

**Ohaus Valor 6000 g x 1g**

Ex demo price £95 AML new price £150



A rapid response IP68 scale with food safety certification (NSF) that also supports an HACCP compliant quality system. It is built to enhance productivity in harsh environments, having large (20.5 mm. high digits) red LED displays front and rear,

Shipping £9.50 per order. Calibration Certificate £15 All prices plus Vat.

Call 01427 668 806 or email [alex@amlinstruments.co.uk](mailto:alex@amlinstruments.co.uk) for more details

**Ex Demo & Used Balance & Scales List.**  
**If you can't find what you need please contact us.**

**Ranger Count 30 kg x 1 g**

Ex demo price £195 AML new price £295



The Ohaus Ranger 3000 counting scale is designed to be the best compact counting scale in its class. It is a high resolution machine able to give stable and accurate readings within one second.

Rugged ABS housing and stainless steel platform

Powered by mains or rechargeable battery (charger included)

Up to 210 hours operation per charge

**UWE APM 150 kg x 0.05 kg Legal for trade approved scale**

Ex demo price £125



The large LCD display with backlight gives excellent readability in sunlight which makes the APM Series ideal for outdoor use. The heavy duty construction will withstand harsh industrial operation yet its light weight makes it ideal for portability and the low profile makes it easy to load and unload the platform.

Shipping £9.50 per order. Calibration Certificate £15 All prices plus Vat.

Call 01427 668 806 or email [alex@amlinstruments.co.uk](mailto:alex@amlinstruments.co.uk) for more details

**Ex Demo & Used Balance & Scales List.**  
**If you can't find what you need please contact us.**

**Mettler PE 360 Delta Range 60g x 0.001 g & 60g to 360g 0.01 g**

Second Hand Price £125



**T scale NHB-600 600g x 0.01g**

Ex Demo Price £95 AML New Price £185



Powered by mains or internal rechargeable battery (50 hrs) - charger supplied  
Backlit LCD display with 16.5 mm high digits  
14 weighing units  
RS-232 interface as standard

Shipping £9.50 per order. Calibration Certificate £15 All prices plus Vat.  
Call 01427 668 806 or email [alex@amlinstruments.co.uk](mailto:alex@amlinstruments.co.uk) for more details

**Ex Demo & Used Balance & Scales List.**  
**If you can't find what you need please contact us.**

**LP Series Floor 300 Kg x 20 g.** Plate Size (mm) : 500 x 600 or 60 kg x 5 g version available

Ex Demo Price £110 AML Price New £165



Extremely economical yet capable general purpose **mains powered** floorstanding digital scales with stainless steel indicator unit and weighplate

**Adam Equipment PW 124 120 g x 0.0001g**

Second Hand Price £250 New Price £650



The PW Precision Analytical Balances are a new generation of scales designed for the most demanding applications in both laboratories and industry alike.

This version of the PW has an internal calibration mass. The scale will automatically re-calibrate it's self for temperature change or on a pre-set time basis.

Shipping £9.50 per order. Calibration Certificate £15 All prices plus Vat.

Call 01427 668 806 or email [alex@amlinstruments.co.uk](mailto:alex@amlinstruments.co.uk) for more details