

Motorised vertical test stand SAUTER TVO



Premium test stand for laboratory applications

Features

- Motorised test stand for tension an compression tests
- Table-top design for comfortable operation
- · Robust design for durable use
- Easy-to-access safety switch-off
- Upper and lower end point, can be set individually
- Automatic or manual operation mode
- Can be used for force gauges up to 500 N (e.g. SAUTER FH-S, not included, for details see page 12)

Technical data

- Maximum tensile and compressive force: 500 N
- Maximum travel length: 300 mm
- Speed accuracy: 2 % of [Max]
- Net weight approx. 25 kg

Accessories

- Digital length measuring device, measuring range 300 mm, readout 0,01 mm, details see page 38, SAUTER LB 300-2.
- Mounting the length measuring device LB onto a SAUTER test stand at the factory, SAUTER LB-A02
- Force-time data transfer software for graphical representation on the PC and data transfer to Microsoft Excel®, SAUTER AFH FAST
- Force-displacement data transfer software with graphic display of the measurement process, only in combination with SAUTER LB, SAUTER AFH FD

STANDARD







Model	Measuring range	Speed range	Max. travelling distance	Dimensions	
SAUTER	[Max] N	mm/min	mm	W×D×H mm	
TVO 500N300.	500	15-300	300	236×428×570	



Pictograms



Adjusting program (CAL):

For quick setting of the instrument's accuracy. External adjusting weight required.



Calibration block:

standard for adjusting or correcting the measuring device.



Peak hold function:

capturing a peak value within a measuring process.



Scan mode:

continuous capture and display of measurements



Push and Pull:

the measuring device can capture tension and compression forces.



Length measurement:

captures the geometric dimensions of a test object or the movement during a test process.



Focus function:

increases the measuring accuracy of a device within a defined measuring range.



Internal memory:

to save measurements in the device memory.



Data interface RS-232:

bidirectional, for connection of printer and PC.



Data interface USB:

To connect the measuring instrument to a printer, PC or other peripheral devices.



Data interface Infrared:

To transfer data from the measuring instrument to a printer, PC or other peripheral devices.



Control outputs (optocoupler, digital I/O):

to connect relays, signal lamps, valves, etc.



Analogue interface:

to connect a suitable peripheral device for analogue processing of the measurements



Statistics

using the saved values, the device calculates statistical data, such as average value, standard deviation etc.



PC Software:

to transfer the measurement data from the device to a PC.



Printer:

a printer can be connected to the device to print out the measurement data.



GLP/ISO record keeping:

of measurement data with date, time and serial number. Only with SAUTER printers



Measuring units:

Weighing units can be switched to e.g. non-metric at the touch of a key. Please refer to website for more details.



Measuring with tolerance range (limit-setting function):

Upper and lower limiting can be programmed individually. The process is supported by an audible or visual signal, see the relevant model



ZERO

ZERO:

Resets the display to "0".



Battery operation:

Ready for battery operation. The battery type is specified for each device.



Rechargeable battery pack:

rechargeable set.



Mains adapter:

230V/50Hz in standard version for EU. On request GB, AUS or USA version available.



Power supply:

Integrated, 230V/50Hz in EU. More standards e.g. GB, AUS or USA on request.



Motorised drive:

The mechanical movement is carried out by a electric motor.



Motorised drive:

The mechanical movement is carried out by a synchronous motor (stepper).



Fast-Move:

the total length of travel can be covered by a single lever movement.



DAkkS calibration possible:

The time required for DAkkS calibration is shown in days in the pictogram.



Factory calibration:

The time required for factory calibration is specified in the pictogram.



Package shipment:

The time required for internal shipping preparations is shown in days in the pictogram.



Pallet shipment:

The time required for internal shipping preparations is shown in days in the pictogram.

Your KERN specialist dealer: